



WOVEN WIRE SIEVES

Endecotts woven wire mesh sieves are the most widely used test sieves for all types of laboratory sampling and particle size analysis.

They are made with only the highest quality materials and are available in diameter sizes of 100, 200, 300, 400 or 450 mm as well as in 8 or 12 inches.

Woven wire mesh sieves can be supplied with aperture sizes ranging from 125 mm down to 20 microns in full and half height versions.

International Test Sieve Series (ISO) American Standard Test Sieve Series (ASTM)



PRODUCT ADVANTAGES

- | Precision frame (ensures consistent nestability)
- | Precise aperture (in accordance with ISO 3310-1, ASTM E11 or other specifications)
- | Natural fillet (free flowing of sample)
- | Totally sealed (no crevice to lose material)
- | Evenly tensioned mesh (ensures accurate analysis)
- | Safe edge (big radius makes it comfortable to handle)
- | Serial number (ensures full traceability)

WOVEN WIRE SIEVES

TECHNICAL DATA

Type of mesh	woven wire
Diameter	100 / 200 / 300 / 400 / 450 mm 8 / 12 inches
Mesh sizes	125 mm - 20 µm
Frame material	stainless steel
Mesh material	stainless steel
Height	full and half (200 / 300 mm 8 / 12 inches) height version



www.endecotts.com/sieves

ORDER DATA

Ø 100 MM STAINLESS STEEL				FULL HEIGHT ISO 3310-1		
	# mm	frame	mesh	standard	height	ø
100SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW14.0	14.00 mm	stainless	stainless	ISO 3310-1	full	100 mm

		steel	steel			
100SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm

		steel	steel			
100SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.900	900 μ m	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.850	850 μ m	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.800	800 μ m	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.710	710 μ m	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.630	630 μ m	stainless steel	stainless steel	ISO 3310-1	full	100 mm

		steel	steel			
100SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm

		steel	steel			
100SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm

		steel	steel			
100SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	100 mm
100S_STLID	Lid	stainless steel	100 mm			
100S_STREC	Receiver	stainless steel	full	100 mm		
100S_STINT	Intermediate Receiver	stainless steel	full	100 mm		

Ø 100 MM STAINLESS STEEL FULL HEIGHT ASTM E11

	# mm	frame	mesh	standard	height	Ø	
100SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	101 mm
100SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	100 mm

100SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	100 mm

100SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	100 mm

100SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	100 mm
100SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	100 mm
100S_STLID	Lid	stainless steel	100 mm				
100S_STREC	Receiver	stainless steel	full	100 mm			
100S_STINT	Intermediate Receiver	stainless steel	full	100 mm			

Ø 200 MM STAINLESS STEEL FULL HEIGHT ISO 3310-1

	# mm	frame	mesh	standard	height	Ø
200SIW125.	125.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW75.0	75.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm

		steel	steel			
200SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm

		steel	steel			
200SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm

		steel	steel			
200SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm

		steel	steel			
200SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm

		steel	steel			
200SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	200 mm
200S_STLID	Lid	stainless steel	200 mm			
200S_STREC	Receiver	stainless	full	200 mm		

steel

200S_STINT	Intermediate Receiver	stainless steel	full	200 mm
------------	-----------------------	-----------------	------	--------

Ø 200 MM STAINLESS STEEL FULL HEIGHT ASTM E11

	# mm	frame	mesh	standard	height	Ø	
200SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	200 mm



200SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	200 mm

200SAW.118	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	200 mm

200SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	200 mm
200SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	200 mm
200S_STLID	Lid	stainless steel	200 mm				
200S_STREC	Receiver	stainless steel	full	200 mm			
200S_STINT	Intermediate Receiver	stainless steel	full	200 mm			

Ø 300 MM STAINLESS STEEL FULL HEIGHT ISO 3310-1

	# mm	frame	mesh	standard	height	Ø
300SIW125.	125.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW75.0	75.00 mm	stainless	stainless	ISO 3310-1	full	300 mm

		steel	steel			
300SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm

		steel	steel			
300SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm

		steel	steel			
300SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm

		steel	steel			
300SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm

		steel	steel			
300SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm

		steel	steel			
300SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	300 mm
300S_STLID	Lid	stainless steel	300 mm			
300S_STREC	Receiver	stainless steel	full	300 mm		
300S_STINT	Intermediate Receiver	stainless steel	full	300 mm		
300S_STRECH	Receiver	stainless steel	half	300 mm		
300S_STINTH	Intermediate Receiver	stainless steel	half	300 mm		

Ø 300 MM STAINLESS STEEL FULL HEIGHT ASTM E11

	# mm	frame	mesh	standard	height	ø	
300SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	300 mm

300SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	300 mm

300SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	300 mm

300SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	300 mm
300SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	300 mm
300S_STLID	Lid	stainless steel	300 mm				
300S_STREC	Receiver	stainless steel	full	300 mm			
300S_STINT	Intermediate Receiver	stainless steel	full	300 mm			

Ø 400
MM

STAINLESS STEEL

FULL
HEIGHT

ISO 3310-1

	# mm	frame	mesh	standard	height	Ø
400SIW125.	125.00 mm	stainless	stainless	ISO 3310-1	full	400 mm

		steel	steel			
400SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW75.0	75.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm

		steel	steel			
400SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm

		steel	steel			
400SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm

		steel	steel			
400SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm

		steel	steel			
400SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm

		steel	steel			
400SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	400 mm
400S_STLID	Lid	stainless steel	400 mm			
400S_STREC	Receiver	stainless steel	full	400 mm		
400S_STINT	Intermediate Receiver	stainless steel	full	400 mm		

Ø 400 MM	STAINLESS STEEL		FULL HEIGHT	ASTM E11			
	# mm	frame	mesh	standard	height	ø	
400SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	400 mm

400SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	400 mm

400SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	400 mm

400SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	400 mm
400SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	400 mm
400S_STLID	Lid	stainless steel	400 mm				
400S_STREC	Receiver	stainless steel	full	400 mm			

400S_STINT Intermediate stainless full 400 mm
Receiver steel

Ø 450 MM STAINLESS STEEL FULL HEIGHT ISO 3310-1

	# mm	frame	mesh	standard	height	Ø
450SIW125.	125.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW75.0	75.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW37.5	37.50 mm	stainless	stainless	ISO 3310-1	full	450 mm

		steel	steel			
450SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm

		steel	steel			
450SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm

		steel	steel			
450SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm

		steel	steel			
450SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm

		steel	steel			
450SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	450 mm
450S_STLID	Lid	stainless steel	450 mm			
450S_STREC	Receiver	stainless steel	full	450 mm		
450S_STINT	Intermediate Receiver	stainless steel	full	450 mm		

Ø 450 MM STAINLESS STEEL FULL HEIGHT ASTM E11

mm frame mesh standard height ø

450SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	450 mm

450SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	450 mm

450SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	450 mm

450SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	450 mm
450SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	450 mm
450S_STLID	Lid	stainless steel	450 mm				
450S_STREC	Receiver	stainless steel	full	450 mm			
450S_STINT	Intermediate Receiver	stainless steel	full	450 mm			

Ø 8" STAINLESS STEEL FULL HEIGHT ISO 3310-1

	# mm	frame	mesh	standard	height	Ø
008SIW125.	125.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW75.0	75.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"

		steel	steel			
008SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	8"

		steel	steel			
008SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	8"

		steel	steel			
008SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	8"

		steel	steel			
008SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	8"

		steel	steel			
008SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	8"
008S_STLID	Lid	stainless steel	8"			
008S_STREC	Receiver	stainless	full	8"		

		steel		
008S_STINT	Intermediate Receiver	stainless steel	full	8"
008S_STRECH	Receiver	stainless steel	half	8"
008S_STINTH	Intermediate Receiver	stainless steel	half	8"

Ø 8" STAINLESS STEEL FULL HEIGHT ASTM E11

	# mm	# mesh number	frame	mesh	standard	height	Ø
008SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	8"
008SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	8"
008SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	8"
008SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	8"
008SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	8"
008SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	8"
008SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	8"
008SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	8"
008SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	8"
008SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	8"
008SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	8"
008SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	8"

			steel	steel			
008SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	8"
008SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	8"
008SAW19.0	19.00 mm	3/4	stainless steel	stainless steel	ASTM E11	full	8"
008SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	8"
008SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	8"
008SAW12.5	12.50 mm	1/2	stainless steel	stainless steel	ASTM E11	full	8"
008SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	8"
008SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	8"
008SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	8"
008SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	8"
008SAW6.30	6.30 mm	1/4	stainless steel	stainless steel	ASTM E11	full	8"
008SAW5.60	5.60 mm	No. 3 1/2	stainless steel	stainless steel	ASTM E11	full	8"
008SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	8"
008SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	8"
008SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	8"
008SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	8"
008SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	8"
008SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	8"

			steel	steel			
008SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	8"
008SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	8"
008SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	8"
008SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	8"

			steel	steel			
008SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	8"
008SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	8"
008S_STLID	Lid	stainless steel	8"				
008S_STREC	Receiver	stainless steel	full	8"			
008S_STINT	Intermediate Receiver	stainless steel	full	8"			
008S_STRECH	Receiver	stainless steel	half	8"			
008S_STINTH	Intermediate Receiver	stainless steel	half	8"			

Ø 12" STAINLESS STEEL FULL HEIGHT ISO 3310-1						
	# mm	frame	mesh	standard	height	Ø
012SIW125.	125.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW112.	112.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW106.	106.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW100.	100.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW90.0	90.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW80.0	80.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW75.0	75.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW71.0	71.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW63.0	63.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW56.0	56.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW53.0	53.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW50.0	50.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW45.0	45.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW40.0	40.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW37.5	37.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW35.5	35.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW31.5	31.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW28.0	28.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW26.5	26.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW25.0	25.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW22.4	22.40 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW20.0	20.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW19.0	19.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW18.0	18.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW16.0	16.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW14.0	14.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW13.2	13.20 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW12.5	12.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW11.2	11.20 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW10.0	10.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW9.50	9.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW9.00	9.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW8.00	8.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW7.10	7.10 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW6.70	6.70 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW6.30	6.30 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW5.60	5.60 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW5.00	5.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW4.75	4.75 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW4.50	4.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW4.00	4.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW3.55	3.55 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW3.35	3.35 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW3.15	3.15 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW2.80	2.80 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW2.50	2.50 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW2.36	2.36 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW2.24	2.24 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW2.00	2.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.80	1.80 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.70	1.70 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.60	1.60 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.40	1.40 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.25	1.25 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.18	1.18 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW1.12	1.12 mm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW1.00	1.00 mm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.900	900 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.850	850 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.800	800 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.710	710 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.630	630 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.600	600 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.560	560 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.500	500 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.450	450 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.425	425 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.400	400 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.355	355 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.315	315 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.300	300 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.280	280 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.250	250 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.224	224 µm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW.212	212 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.200	200 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.180	180 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.160	160 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.150	150 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.140	140 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.125	125 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.112	112 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.106	106 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.100	100 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.090	90 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.080	80 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.075	75 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.071	71 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.063	63 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.056	56 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.053	53 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.050	50 µm	stainless steel	stainless steel	ISO 3310-1	full	12"

012SIW.045	45 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.040	40 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.038	38 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.036	36 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.032	32 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.025	25 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012SIW.020	20 µm	stainless steel	stainless steel	ISO 3310-1	full	12"
012S_STLID	Lid	stainless steel	12"			
012S_STREC	Receiver	stainless steel	full	12"		
012S_STINT	Intermediate Receiver	stainless steel	full	12"		
012S_STRECH	Receiver	stainless steel	half	12"		
012S_STINTH	Intermediate Receiver	stainless steel	half	12"		

Ø 12" STAINLESS STEEL FULL HEIGHT ASTM E11							
	# mm	# mesh number	frame	mesh	standard	height	Ø
012SAW125.	125.00 mm	5.00	stainless steel	stainless steel	ASTM E11	full	12"
012SAW106.	106.00 mm	4.24	stainless steel	stainless steel	ASTM E11	full	12"
012SAW100.	100.00 mm	4	stainless steel	stainless steel	ASTM E11	full	12"

012SAW90.0	90.00 mm	3 ½	stainless steel	stainless steel	ASTM E11	full	12"
012SAW75.0	75.00 mm	3	stainless steel	stainless steel	ASTM E11	full	12"
012SAW63.0	63.00 mm	2 ½	stainless steel	stainless steel	ASTM E11	full	12"
012SAW53.0	53.00 mm	2.12	stainless steel	stainless steel	ASTM E11	full	12"
012SAW50.0	50.00 mm	2	stainless steel	stainless steel	ASTM E11	full	12"
012SAW45.0	45.00 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	12"
012SAW37.5	37.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	12"
012SAW31.5	31.50 mm	1 ¾	stainless steel	stainless steel	ASTM E11	full	12"
012SAW26.5	26.50 mm	1.06	stainless steel	stainless steel	ASTM E11	full	12"
012SAW25.0	25.00 mm	1	stainless steel	stainless steel	ASTM E11	full	12"
012SAW22.4	22.40 mm	7/8	stainless steel	stainless steel	ASTM E11	full	12"
012SAW19.0	19.00 mm	¾	stainless steel	stainless steel	ASTM E11	full	12"
012SAW16.0	16.00 mm	5/8	stainless steel	stainless steel	ASTM E11	full	12"
012SAW13.2	13.20 mm	0.530	stainless steel	stainless steel	ASTM E11	full	12"
012SAW12.5	12.50 mm	½	stainless steel	stainless steel	ASTM E11	full	12"
012SAW11.2	11.20 mm	7/16	stainless steel	stainless steel	ASTM E11	full	12"
012SAW9.50	9.50 mm	3/8	stainless steel	stainless steel	ASTM E11	full	12"
012SAW8.00	8.00 mm	5/16	stainless steel	stainless steel	ASTM E11	full	12"

012SAW6.70	6.70 mm	0.265	stainless steel	stainless steel	ASTM E11	full	12"
012SAW6.30	6.30 mm	¼	stainless steel	stainless steel	ASTM E11	full	12"
012SAW5.60	5.60 mm	No. 3 ½	stainless steel	stainless steel	ASTM E11	full	12"
012SAW4.75	4.75 mm	No. 4	stainless steel	stainless steel	ASTM E11	full	12"
012SAW4.00	4.00 mm	No. 5	stainless steel	stainless steel	ASTM E11	full	12"
012SAW3.35	3.35 mm	No. 6	stainless steel	stainless steel	ASTM E11	full	12"
012SAW2.80	2.80 mm	No. 7	stainless steel	stainless steel	ASTM E11	full	12"
012SAW2.36	2.36 mm	No. 8	stainless steel	stainless steel	ASTM E11	full	12"
012SAW2.00	2.00 mm	No. 10	stainless steel	stainless steel	ASTM E11	full	12"
012SAW1.70	1.70 mm	No. 12	stainless steel	stainless steel	ASTM E11	full	12"
012SAW1.40	1.40 mm	No. 14	stainless steel	stainless steel	ASTM E11	full	12"
012SAW1.18	1.18 mm	No. 16	stainless steel	stainless steel	ASTM E11	full	12"
012SAW1.00	1.00 mm	No. 18	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.850	850 µm	No. 20	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.710	710 µm	No. 25	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.600	600 µm	No. 30	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.500	500 µm	No. 35	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.425	425 µm	No. 40	stainless steel	stainless steel	ASTM E11	full	12"

012SAW.355	355 µm	No. 45	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.300	300 µm	No. 50	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.250	250 µm	No. 60	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.212	212 µm	No. 70	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.180	180 µm	No. 80	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.150	150 µm	No. 100	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.125	125 µm	No. 120	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.106	106 µm	No. 140	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.090	90 µm	No. 170	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.075	75 µm	No. 200	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.063	63 µm	No. 230	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.053	53 µm	No. 270	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.045	45 µm	No. 325	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.038	38 µm	No. 400	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.032	32 µm	No. 450	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.025	25 µm	No. 500	stainless steel	stainless steel	ASTM E11	full	12"
012SAW.020	20 µm	No. 635	stainless steel	stainless steel	ASTM E11	full	12"
012S_STLID	Lid	stainless steel	12"				

012S_STREC	Receiver	stainless steel	full	12"
012S_STINT	Intermediate Receiver	stainless steel	full	12"
012S_STRECH	Receiver	stainless steel	half	12"
012S_STINTH	Intermediate Receiver	stainless steel	half	12"