



WET WASHING SIEVES

Wet washing sieves are extremely useful where samples need to be separated with the help of liquids. These sieves are available in 200 mm or 8 inch diameter (4 or 8 inches deep), with brass or stainless steel frames.

Endecotts supplies wet washing sieves for a wide range of aperture sizes with optional support medium for fine mesh.



PRODUCT ADVANTAGES

- | Precision frame (ensures consistent nestability)
- | Precise aperture (in accordance with ISO 3310, ASTM or other specifications)
- | Natural fillet (free flowing of sample)
- | Totally sealed (no crevice to lose material)
- | Evenly tensioned mesh (ensures accurate analysis)
- | Safe edge (big radius makes it comfortable to handle)
- | Serial number (ensures full traceability)

WET WASHING SIEVES

TECHNICAL DATA

Type of mesh	woven wire
Diameter	200 / 450 mm / 8 inches
Mesh sizes	125 mm - 20 μ m
Frame material	stainless steel
Mesh material	stainless steel
Height	100 / 200 mm 4 / 8 inches 300mm for 450 Dia.

www.endecotts.com/sieves



ORDER DATA